

DSA 1000 unit

Technical data



Voltage: 230 V 50/60 Hz

Power rating: 350 W

Main fuses: 2x8 A

Compressor: 15 l/min. 3 bar

Reservoir capacity: 2 x 1 l

Dimensions HxWxD (cm): 79x39x55

Colour: RAL 9010

Corian instrument tray available in 9 colours.

Instruments from left to right (front view):

- * Position 1: EMS scaler or BIEN AIR MC40 motor
- * Position 2: EMS scaler or BIEN AIR MC3 motor
- * Position 3: LUZZANI multi-function syringe