

DSA 2000 L unit

Technical data



Voltage: 230 V 50/60 Hz

Power rating: 400 W

Main fuses: 2x8 A

Compressor: 50 l/min. 8 bar

Reservoir capacity: 2 x 2 l

Dimensions HxWxD (cm): 79x45x62

Colour: RAL 9002

Corian instrument tray available in 9 colours.

Instruments from left to right (front view):

- * Position 1: FARO multi-function syringe
- * Position 2: cord for turbine
- * Position 3: BIEN AIR MC40 micro motor
- * Position 4: BIEN AIR MC40 micro motor
- * Position 5: EMS or SATELEC scaler