

DSA 2500 unit

Technical data



Voltage: 230 V 50/60 Hz

Power rating: 450 W

Main fuses: 2x8 A

Compressor: 50 l/min. 8 bar

Reservoir capacity: 2 x 2 l

Dimensions HxWxD (cm): 79x45x60

Colour: RAL 9002

Corian instrument

Instruments from left to right (front view):

- * Position 1: FARO multi-function syringe
- * Position 2: cord for turbine or LUZZANI syringe or BIEN AIR MC3
- * Position 3: BIEN AIR MC3 micro motor or cord for turbine
- * Position 4: BIEN AIR MC3 micro motor or cord for turbine
- * Position 5: cord for turbine or BIEN AIR MC3 or LED curing light
- * Position 6: EMS or SATELEC scaler